

dream²





Sensors monitor your sleep with lab-level precision.

Active audio features improve your sleep all night long.

Daily reports use the data collected by the headband for a clear, complete vision of your night.

Coaching programs designed to help you resolve your sleep problems.

The Dreem solution




The headband

- Lab grad sleep analysis
- Biofeedback audio relaxation techniques
- Deep sleep stimulations
- Smart Alarm



The mobile app

- Daily sleep reports
- Tracking of the progress
- Sleep restructuring program



Dreem is a complete, personalized non-invasive and non-chemical solution for people suffering from sleep disorders and insomnia.



dreem²

No Wi-fi, no Bluetooth
during the night

Supple materials and foam
for more flexibility and
comfort.

Thanks to its set of
adjusters (sizes S,M
and L), Dreem 2
adapts to all different
head sizes.




A redesigned
touchpad, that is
easier to navigate
during the night.

Sensors that monitor
your sleep with
precision.

Active audio features,
transmitted by bone
conduction.

A sleep clinic at home



Dreem monitors sleep with the same accuracy as a polysomnography, the tool used in sleep clinic.



* <https://www.biorxiv.org/content/10.1101/662734v1>

The Dreem app enables you to **understand and improve your sleep.**

Every morning, **you will have access to your sleep data** and the app will provide you with **weekly reports** to measure your progress.



The Dreem sleep restructuring program

Dreem has developed a tailored program based on the reference treatment for insomnia which is more engaging and efficient than a classic therapy in sleep clinics.

A more **engaging** and **efficient** program based on accurate and objective data.

Developed with **sleep experts** with a daily support through the app **directly from home.**

Dreem supports you during your entire program and you can **book calls with our sleep experts.**



Biofeedback audio relaxation techniques



Ambiance



Meditation



Cognition



Respiration



Pink noises



Smart alarm

- **Falling asleep**

Dreem proposes active audio features to fall asleep faster.

- **Deep sleep**

Deep sleep stimulations, with scientifically proven efficacy will enhance the quality of your deep sleep.

- **Waking up**

Dreem wakes you up at the optimal moment to start the day refreshed.

Falling asleep

Dreem offers a series of audio relaxation techniques to help you **fall asleep more easily.**



Ambiance



Meditation



Cognition



Breathing



Deep sleep

Sound stimulation, of which the effectiveness is scientifically proven, **improve the quality of your deep**



Pink Noise



Falling asleep


Dreem wakes you up at the optimal time for an easier morning and less tiredness during the day.




Smart alarm




Renowned sleep specialists support Dreem.




Chief Scientist Officer & President of the Allen Institute for Brain Sciences, Seattle, US.



Stanford Sleep Center Director



Co-Director of CIRCS, Lausanne, Switzerland.



Head of the Sleep and Circadian Neuroscience Institute



Our **users** say it best:

dreem

Dreem

618 Reviews - Good



I've tried lots of trackers, but this is the only one that I feel is giving me an accurate measure of my sleep because with other ones, I can just lay there wide awake and they think I am sleeping. I was scared the headband wouldn't be comfortable, but it hasn't bothered me, **the EEG makes it feel like I am at a sleep clinic in my own bed.** I just finished my assessment a few days ago and have been following their coaching. So far I've been following along and it seems like **I am making some improvements based on how I feel in the morning and what the data shows.** Highly recommend to anyone who is struggling with insomnia or any issues with their sleep.



Jaybee

Dreemer for 1 year

In short, the Dreem headband has been helping. The effect is slight but noticeable, and I awaken each morning for the first time in years with some small feeling of restoration. Given that I **have long complained of suffering from non-restorative sleep, this device has been a godsend.**



Foe

Dreemer for 7 months

Like me before Dreem, you can't sleep good ? Your health is impacted by a low sleep quality ? Your life quality is impacted too ? **You have to try Dreem and coach your brain to learn how to find a good quality of sleep. It's not a miracle : Dreem and time are your pals for a better sleep.** It's only a revolution !



Pierre

Dreemer for 9 months

Better sleep made
even more accessible.

Dreem 2
399€ (RRP)



Weight 130g	Size Adjustable with S, M, L band size adjusters	Dimension Head circumference 520mm - 610mm	Interface Capacitive sensor RGB LED indicator
Battery 1200mAh Up to 12 hours	Charging time 2h30	Sensors EEG sensors x 6 Pulse oximeter Accelerometer Sonometer	Materials Hypoallergenic TPU Soft fabric
Audio Bone Conductor system 3.5mm Audio Jack	Connectivity Bluetooth 4.0 Low Energy Bluetooth Classic Micro USB Wifi b/g/n Pogo charging pins	Compatibility Android 5.0 and more recent iOS 10 and more recent iPhone 5 and more recent iPod Touch 6th Gen and more recent	Languages English French

dreem²